
CareVis®

EINMALINFUSIONSPUMPEN / SINGLE USE INFUSION PUMP

GER

ENG

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SINGLE USE INFUSION PUMP

Description

CareVis® is a single-use, portable infusion pump for patients which works without batteries or electricity. It provides liquid medication at an accurate continuous pre-fixed flow with either one (CareVis®) or multiple (CareVis® Multirate) selectable rates or, if selected, additional bolus doses via a patient-controlled analgesia (PCA) module (CareVis® PCA). It consists of an elastomeric balloon inside a rigid, transparent container and an infusion line with one or two (CareVis® Twin) flow restrictors and Luer-Lock patient connectors. The PCA module has a bolus button. When pressed, medication in a bag beneath the bolus button is transferred to the patient. The filling time (= lock-out time) of this bag is pre-fixed by a secondary flow restrictor. The Multirate module has a detachable key, which allows switching between different continuous flow rates. The protective transparent case of the device provides UV blocking capabilities up to 390nm, blocking UVB, UVC and most UVA light. CareVis® is equipped with a particle filter (12) of 1.2µm and a 0.03 µm pore-size air filter which forces air bubbles out of the system.

Indications for use

CareVis® is intended for the infusion of chemotherapy, antibiotic therapy, analgesic therapy and other general intravenous, intraarterial, subcutaneous, epidural, perineural and intraligamentary (in operative wounds) infusions. CareVis® PCA, and CareVis® Multirate have been mainly designed for systemic or regional pain management infusion use. General pain management infusion includes preoperative and postoperative surgery or chronic pain. All CareVis® models are suitable and safe for use with children. However, the pumps should only be used for children and adult patients who, following instruction by medical personnel, fully understand and comprehend the function of the pump and, if chosen, of the bolus button or Multirate unit. If CareVis® Multirate is used, the detachable continuous flow-rate controller may be left on the Multirate (13d) module or removed after setting the desired flow rate. Where the flow-rate controller is left on the module, the patient must be instructed on its use by the medical personnel and must understand it.

Contraindications

Solutions with high viscosity or solid particles should not be used with CareVis®, as they flow slowly and may block the system. As a result, blood, blood products, proteins and fluids containing lipids cannot be used.

Warnings

1. CareVis® is a single-use product: DO NOT REUSE and DO NOT re-sterilise. Reuse or re-sterilisation may result in unpredictable infusion flow and may cause patient infections.

2. After filling CareVis®, reconfirm the information and the flow-rate setting prior to infusion, as incorrect flow rate or other parameters may fatally injure the patient.

3. After filling CareVis® PCA, remove the red holder (19) from the PCA module before use. If the red holder is not removed, flow through the bolus reservoir will not be held back but will be added to the continuous flow. The resulting high volume of infusion medication may cause serious and even fatal injury to the patient.

4. Epidural or perineural administrations of drugs other than those indicated for epidural or perineural use may result in serious irreversible injury and pain to the patient.

5. Medical staff must ensure that a suitable alternative device for the infusion of critical medicines is used when stopping or interrupting the infusion, or changing the infusion rate, might result in serious health problems and even death.

Precautions

CareVis® will work correctly as long as the following precautions are observed.

1. Read the Instructions for Use carefully and don't use CareVis® if you have not understood all the details.

2. CareVis® should be protected from sunlight and UV rays. Keep dry and store at between 0° and 30 °C.

3. Do not use the CareVis® if the packaging is open or defective. Do not re-sterilise.

4. An aseptic technique must be used throughout the filling procedure.

5. The capillary element (13a), the PCA module (13c) and the Multirate module (13d) may deteriorate if placed in contact with solvents. Use water to clean these items.

6. Make sure the infusion line is primed before use. Otherwise, air may be introduced into the patient.

7. The air-filter hole (11) on the air and particle filter (12) must not be covered and must be kept dry and clean. Solvents such as alcohol may damage its seal against the infusion fluid.

8. Do not set the Multirate module (13d) regulator (17) at intermediate positions (between marks) as infusion will stop.

9. During infusion, visually check the progress indicator (5) while the container (4) gradually empties and ensure that infusion is duly proceeding.

10. The pump progress scale on the protective case does not represent the exact volume in the elastomeric reservoir in ml. It is only intended to follow the pump progress by the patient or its caregiver.

11. No elements should be introduced through the air-inlet hole (2) of the container (4).

12. CareVis® can come into contact with water, but no liquid must enter the container (4) through the air-inlet hole (2).

13. Infusion should be interrupted by clamping (10) the line if the pump body or its elements break or leak. Check the capillary element (13a), the PCA module (13c), the Multirate module (13d), the air and particle filter (12) or any section of the infusion line. If CareVis® Twin with 2 outlets is used, both of outlets must be clamped.

14. A small amount of solution will be left in the CareVis® after the infusion has finished (see Table). Any liquid remaining in the infusion line or the system should not be reused.

15. At the end of infusion, dispose of the device according to the local standards on potentially toxic products.

16. When delivering through a catheter system, refer to the instructions provided by the manufacturer for proper use of the system components. The length, inner diameter, and position of the catheter during use can affect nominal flow rate. The use of catheters with a small inner diameter (i.e. 23 Gauge or narrower) is discouraged.

Calibration conditions and information on correct use

1. CareVis® is calibrated under the following conditions:

- The reservoir is filled at nominal volume + residual volume, so that the nominal volume is infused at nominal flow rate.
- The viscosity of the infusion solution used is the same as normal saline (0.9% NaCl solution).
- The temperature of the solution in the container is at room temperature (22°C).
- The temperature of the flow restrictor (capillary element (13a), the capillary tube (13b), the PCA device (13c) or the Multirate module (13d)) is 32 °C, equivalent to the temperature when in contact with the skin.
- The capillary tube (13b) is straight (i.e. not twisted).
- The container (4) and the patient connection (14) are at the same level, with free output.

2. A variation in viscosity has an indirectly proportional effect on the nominal flow rate; an increase/decrease in viscosity slows down/increases the nominal flow rate.

3. Keeping the restrictor in direct contact with the patient's skin helps to stabilise the operating temperature. A variation in temperature has a directly proportional effect on the actual flow rate. Flow rate will decrease/increase approximately 2% per 1°C (1.8°F) decrease/increase in temperature.

4. Lifting the reservoir above the distal connector increases the pressure gradient and the actual flow rate; holding it below has the opposite effect.

Preparing the medication

1. The medication to be infused must be prescribed by the patient's doctor.
2. To prepare the solution, follow the specific instructions for the medication in question, using an aseptic technique.
3. Choose the most suitable CareVis® model for the volume and the flow rate as indicated on the identification label (3). This applies also to bolus volume and lock-out time if a PCA model is selected and to the right Multirate model.

Filling instructions (read carefully before filling)

1. An aseptic technique must be used throughout the procedure.
2. Do not remove the CareVis® from its packaging until ready for use. Once the packaging has been opened, check that all components are in perfect condition and that the infusion line is properly closed with the blue infusion-line end cap (15).
3. While holding the CareVis®, remove the container-filling port cap (9) and retain it for later use.
4. To determine the volume of solution to be infused, it is important to factor in the residual volume (as indicated in Volume Table). Calculate the fill volume by multiplying the desired infusion time (hours) by the nominal flow rate (ml/h). Then, add the residual volume; e.g. to prepare medication with a flow rate of 10ml/hr for 24hours, the filling volume should be 245ml (10ml/hr x 24hr + 5ml).
5. Fill a Luer-Lock syringe with the solution. Insert the syringe into the filling port (8) by turning it gently and locking it to the port. A needle or other sharp instrument should never be used. Fill the elastomeric balloon (6), if needed, by repeating the operation.

During filling, medication flows throughout the tubing, automatically expelling the air from the line through the in-line air filter (12) and the blue end cap filter (15), which has a hydrophobic membrane preventing medication from flowing out. The time required to prime the line should not exceed the time indicated in Priming Time Table.

6. The filling procedure may vary if cytotoxic drugs or drugs with a tendency to precipitation such as 5-Fluorouracil are used. With cytotoxic drugs, it is generally recommended that the patient line be primed first with normal saline and then the pump be filled with the active ingredient immediately afterwards. Fill some normal saline first. When the solution has reached the patient connector (14) – this may take a while, as indicated in the Filling Time Table - clamp the line to prevent osmotic mixture at its distal part. Initial priming of the line with normal saline prevents tissue damage in the event of initial opening extravasation also impedes precipitates within and slowing down of the restrictor module.
7. Partial filling is allowed but may cause a higher flow rate than the nominal rate, i.e. approx. 8% with a reservoir filled at 70% of nominal volume.
8. The maximum volume of the CareVis® indicated in the Volume Table should not be exceeded.
9. The solution should be introduced at a constant speed, without sudden spurts. More force is required when the balloon (6) first starts filling.
10. Check visually that the balloon (6) expands symmetrically along the inner guide of the container (4). If it does not, the device may be defective and should be discarded.
11. Disconnect the syringe and close the filling port (8) firmly with the original cap (9).
12. Complete a patient label (one is included with the product) with all the necessary information and attach it to the container (4).
13. Clamp the line and insert the CareVis® into the bag included with the product for transport, in order to protect the device.

14. IF YOU USE THE CareVis® PCA:

- While filling, do not remove the red holder (19) from the PCA module and do not press the bolus button. Otherwise, air inside the bolus bag cannot be entirely expelled, which may cause incomplete priming of the bolus bag and discharge of residual air.
- After filling, remove the red holder from the PCA module. Don't connect the device to the patient if the holder is not removed. See warnings.

15 IF YOU USE THE CareVis® MULTIRATE

Before filling, make sure that the regulator (17) is set to the maximum flow rate position to accelerate priming of the line. At the end of the filling procedure, set the starting flow rate for the patient or to zero flow rate for later setting.

16. IF YOU USE THE CareVis® TWIN

Before clamping, check that both end lines are completely filled.

Application to patient

1. Before starting infusion, check that the filled device is at room temperature. If the unit has been refrigerated, it should be returned to room temperature before use, which may take up to 6 hours. Do not microwave. Dew may form on the inner surface of the case during this time but it does not adversely affect normal use of the unit. Dew may also form with minor changes of temperature.

2. Before starting infusion, remove the blue end cap from the distal connector, open the clamp and visually verify the appearance of solution at the tip of distal connector.

3. Close the clamp again and connect the distal connector firmly to the device used for the infusion to patient, open the clamp and start infusion.

4. If the Multirate module (13d) is used, turn the regulator (17) to the desired flow-rate position, using the key (18) supplied. Once this operation is complete, remove the key from the regulator and retain it for future use. Important: If the regulator (17) is set to an intermediate position (between marks), infusion will stop.

5. Attach the flow restrictor (capillary element (13a), capillary tube (13b), PCA device (13c), or Multirate device (13d)) to the patient's skin. Ensure that the air and particle filter (12) is kept dry and the air-filter hole (11) is not obstructed.

6. Check the patient label.

7. The infusion is finished when the level indicator (5) is close to 0 and the balloon (6) is empty and fully deflated.

8. Visually inspect the device once the infusion has finished. Clamp (10) the infusion line and disconnect from the patient. To avoid possible handling risks, connect the end of the output connection (14) to the Luer-Lock filling port (8).

9. Dispose of the device according to the centre protocol or local standards.

Information for patients

1. The scale on the pump case is used to dynamically monitor the medication infusion process. If this does not change throughout the infusion, the healthcare provider should be informed about a problem with the flow.

2. Avoid bringing products such as soap or alcohol in contact with the filter and capillary elements, as this may cause leakage from the air-eliminating vent or seals.

3. Do not disconnect the device before infusion is complete unless directed otherwise by the healthcare provider. The infusion is finished when the level indicator (5) is close to 0 and the balloon (6) is empty and fully deflated (tightly attached to the underlying piston).

4. Where the CareVis® PCA is used, the PCA bolus module provides an extra dose of pain medication if required. The bolus button should always be pressed completely flat to the ground without skew and then released.

5. The lock-out time of the CareVis® PCA is the required time to fill the bolus bag, which is located below the bolus button. Pressing the bolus button before the designated lock-out time will cause a reduced bolus dose as the bolus bag will not be completely filled. Pressing the bolus button repeatedly within the lock-out time will not deliver more medication than the designated total bolus volume of medication.